



OPADRY[®] ENTERIC
FORMULATED FOR ORGANIC SOLVENT COATING

Product
Information

Opadry[®] Enteric Product Information

Opadry Enteric is a family of fully formulated, delayed release coating systems for solid oral dosage forms, which are applied by organic or hydro-alcoholic processing techniques.

Specific Opadry Enteric coating formulations have been developed from a choice of enteric polymers, with solubility as a function of the environmental pH in the gastro-intestinal tract. The various Opadry Enteric products available from Colorcon are listed below along with the constituent enteric polymer and pH at which the system starts to dissolve.

Opadry Enteric	Polymer	Polymer Dissolution pH
▪ 91 series	Polyvinyl Acetate Phthalate (PVAP)	pH 5.0
▪ 94 series	Methacrylic Acid Copolymer (Methacrylic acid- methyl methacrylate 1:1 copolymer)	pH 6.0
▪ 95 series	Methacrylic Acid Copolymer (Methacrylic acid- methyl methacrylate 1:2 copolymer)	pH 7.0

OPADRY ENTERIC - 91 SERIES

Opadry Enteric- 91series is formulated using the delayed-release polymer Polyvinyl Acetate Phthalate NF (PVAP, Phthalavin).

Features and Benefits

Features	Benefits
<ul style="list-style-type: none">▪ Fully formulated system▪ Less susceptible to hydrolysis than other phthalate esters such as CAP or HPMCP¹▪ High coating solids concentration (15% coating solids using a hydro-alcoholic solvent system; 5-8% coating solids for an organic solvent system)▪ Pigmented or clear systems are available▪ Simple to dispense and disperse in solvent systems (45 minutes at low shear)▪ Enteric protection is achieved with 6 - 8% coating weight gain on tablets**▪ Polymer start to dissolve at pH 5.0 and above	<ul style="list-style-type: none">▪ Time and cost savings▪ Stable product (12 months re-evaluation period*)▪ Produces elegant glossy coatings while exhibiting excellent enteric protection▪ Allows rapid coating application▪ Provides rapid, reproducible drug release in alkaline dissolution media.

Regulatory

- Acceptable for use in pharmaceutical enteric coatings in the U.S. and other countries that accept NF compliance. Not acceptable for use in dietary supplement coatings in the U.S.
- Currently no precedence of use for PVAP in Japan exists for pharmaceutical or dietary supplement coatings.

Applications

- Tablets, soft-gelatin capsules, hard-gelatin capsules and multi-particulates.

Packaging

- 25 kg carton with double polyethylene bags and desiccant between bags.

Recommended Storage Conditions

- Seal container tightly.
- Keep below 30°C / 86°F; not more than 65% relative humidity. Product packaging contains desiccant.

Re-evaluation Period (Warranty)

- 12- months from date of manufacture*

* May vary based on pigment selection.

** To be determined on a case-by-case basis.

OPADRY ENTERIC - 94 SERIES

Opadry Enteric- 94series is formulated using the delayed-release polymer

“Methacrylic Acid Copolymer, Type A” USP/ NF

“Methacrylic Acid- Methyl Methacrylate Copolymer (1:1)” Ph. Eur

“Methacrylic Acid Copolymer L” JPE

Features and Benefits

Features	Benefits
<ul style="list-style-type: none">▪ Fully formulated system▪ Non-tacky▪ Hydro-alcoholic solvents are recommended▪ Pigmented or clear systems are available▪ Simple to dispense and disperse in solvent systems (45 minutes at low shear)▪ Enteric protection is achieved with 5 - 6% coating weight gain on tablets**▪ Polymer start to dissolve at pH 6.0 and above▪ Methacrylic acid copolymers have excellent stability against hydrolysis²	<ul style="list-style-type: none">▪ Time and cost savings▪ Processing ease▪ Allows rapid coating application▪ Produces elegant glossy coatings while exhibiting excellent enteric protection▪ Provides rapid, reproducible drug release in alkaline dissolution media.▪ Stable product (12 months re-evaluation period*)

Regulatory

- Acceptable for use in pharmaceutical enteric coatings in the U.S/ Europe and Japan. Not acceptable for use in dietary supplement coatings in the U.S.

Applications

- Tablets, soft-gelatin capsules, hard-gelatin capsules and multi-particulates.

Packaging

- 25 kg carton with double polyethylene bags and desiccant between bags.

Recommended Storage Conditions

- Seal container tightly.
- Keep below 30°C / 86°F; not more than 65% relative humidity. Product packaging contains desiccant.

Re-evaluation Period (Warranty)

- 12- months from date of manufacture*

* May vary based on pigment selection.

** To be determined on a case-by-case basis.

References

1. Hogan, J. E. (1995). Modified Release Coatings. In: G. Cole (Ed.) *Pharmaceutical Coating Technology* (p. 435). London: Taylor & Francis Ltd.
2. Lehmann, K. (1997). Polymethacrylate Coating Systems. In: McGinity J (Ed.) *Aqueous Polymeric Coatings for Pharmaceutical Dosage Forms* (p.136) New York: Marcel Dekker Inc.

World Headquarters

Colorcon

415 Moyer Blvd., P.O. Box 24, West Point, PA 19486-0024

Tel: 215-699-7733 Fax: 215-661-2605 Website: www.colorcon.com

Locations	Telephone	Facsimile	Locations	Telephone	Facsimile
<i>United States</i>			<i>Asia/Pacific</i>		
Irvine, California	949-753-5780	949-450-0452	Singapore	65-6438-0318	65-6438-0178
Indianapolis, Indiana	317-545-6211	317-545-6218	Fuji-gun, Shizuoka, Japan	81-5-4465-2711	81-5-4465-2730
Humacao, Puerto Rico	787-852-3815	787-852-0030	Shanghai, China	86-21-5442-2222	86-21-5442-2229
Stoughton, Wisconsin	608-887-8970	608-887-8984	Goa, India	91-832-288-3434	91-832-288-3440
<i>Canada</i>			Gyeonggi-do, Korea	82-31-296-2173	82-31-296-2178
St. Laurent, QC, Canada	514-337-8341	514-337-9159	<i>Latin America</i>		
<i>Europe</i>			Buenos Aires, Argentina	54-11-4552-1565	54-11-4552-3997
Dartford, Kent, England	44-1322-293000	44-1322-627200	Cotia, Brasil	55-11-4612-4262	55-11-4612-3307
Massy, France	33-1-6447-9750	33-1-6932-5983	Bogota, Colombia	571-418-1202	571-418-1257
Idstein, Germany	49-6126-9961-0	49-6126-9961-11	Caracas, Venezuela	58-212-237-9842	58-212-238-2259
Gallarate, Italy	39-0331-776932	39-0331-776831	Mexico DF, Mexico	52-55-3000-5700	52-55-3000-5701 /02
Budapest, Hungary	36-1-200-8000	36-1-200-8010			
Istanbul, Turkey	90-216-465-0360	90-216-465-0361			
Barcelona, Spain	34-9-3589-3756	34-9-3589-3792			

The information contained herein, to the best of our knowledge is true and accurate. Any recommendations or suggestions are made without warranty or guarantee, since the conditions of use are beyond our control. Any information contained herein is intended as a recommendation for use of our products so as not to infringe on any patent.

Copyright © 2008 Colorcon. The information contained in this document is proprietary to Colorcon and may not be used or disseminated inappropriately.

All trademarks, except where noted, are property of BPSI Holdings, LLC.
mr/opadry_enteric/prod.info/01.2009