

Complete Film Coating System for Nutritionals

Coating Process Conditions – Clear Formula

Opadry[®] NS is a complete film coating system for nutritional supplements. The coating parameters listed in the table below are based on Colorcon data. Individual product and machine functions should be taken into account and conditions altered as required.

	Fully	Perforated Side Ven	ted Coating Pan		
Pan Diameter	cm	38.1	60.96	122	152.4
Fall Diameter	in	15	24	48	60
		Solution Prepa	ration		
Solvent		Purified Water	Purified Water	Purified Water	Purified Wat
Solids Content	% w/w	10	10	10	10
		Run Objectiv			
Theoretical Weight Gain	% w/w	1	1	1	1
Tablet Charge	kg	1.5 - 3.0	12 - 15	120 - 160	320 - 350
		Gun Configura			
Example Spray Gun	Туре	Schlick 931/7-1 S22	Spraying Systems 1/8 VAU	Schlick 931/7-1 S35 (ABC)	Manesty
Number of Guns	guns	1	1-2	3 - 4	5 - 6
Fluid Nozzle	type	Standard	VF4078-316SS	Standard	Standard
Nozzle Diameter	mm	1.2	1.02 mm	1.2	1.2
	in	0.05	0.04 in	0.05	0.05
Air Cap	type	Standard (#46 - Horned)	VA113293-60- 316SS	ABC - Anti Bearding Cap	Standard (Horned)
Gun to Bed Distance	cm	10	17	25	25
	in	4	7	10	10
Gun to Gun Distance	cm	N/A	10	15	15
	in	N/A	4	6	6
Atomizing Air Pressure	bar	1.4	1.7	2.5	2.5
	psi	20	25	36	36
Pattern Air Pressure	bar	1.4	1.7	2.5	2.5
	psi	20	25	36	36
		Process Param	neters		
Drying Air Volume	cu.m/hr	250 - 290	424 - 470	2300 - 3060	5440 - 663
	cu.ft/min	147 - 175	250 - 275	1350 - 1800	3200 - 390
Pan Speed	rpm*	18	14	8	6
Pre-warm tablet Bed	°C	45 - 50	45 - 50	45 - 50	45 - 50
	°F	113 - 122	113 - 122	113 - 122	113 - 122
Inlet Air Temperature	°C	60	65	65	65
	°F	140	149	149	149
Tablet Bed Temperature	°C	43 - 48	43 - 48	43 - 48	43 - 48
	°F	109 - 118	109 - 118	109 - 118	109 - 118
Spray Rate	g/min	12 - 20	40 - 60	260 - 430	650 - 750
Approximate Coating Times	min	15 - 20	30 - 40	45 - 60	55 - 65

*Pan speed should be set so that even mixing of the bed occurs.

This will be influenced by tablet size, shape, friability and baffle type in the pan.



The information contained herein, to the best of Colorcon, Inc.'s knowledge is true and accurate. Any recommendations or suggestions of Colorcon, Inc. with regard to the products provided by Colorcon, Inc. are made without warranty, either implied or expressed, because of the variations in methods, conditions and equipment which may be used in commercially processing the products, and no such warranties are made for the suitability of the products for any applications that you may have disclosed. Colorcon, Inc. shall not be liable for loss of profit or for incidental, special or consequential loss or damages.

Colorcon, Inc. makes no warranty, either expressed or implied, that the use of the products provided by Colorcon, Inc., will not infringe any trademark, trade name, copyright, patent or other rights held by any third person or entity when used in the customer's application.

For more information, contact your Colorcon representative or call:

North America +1-215-699-7733 Europe/Middle East/Africa +44-(0)-1322-293000 Asia Pacific +65-6438-0318 Latin America +54-11-4552-1565



© BPSI Holdings LLC, 2010. The information contained in this document is proprietary to Colorcon, Inc. and may not be used or disseminated inappropriately.

All trademarks, except where noted, are property of BPSI Holdings, LLC.