

Aqueous Acrylic Enteric System

# Optimal Coating Process Parameters for a, Fully Formulated, Acrylic-Based, Enteric, Film Coating System

# OBJECTIVES

- Identify the sensitivity of each response variable to the independent process variables.
- Determine the optimal coating process parameters.
- Demonstrate the overall applicability of a developmental film coating system.

# METHODOLOGY

### Materials

The developmental film coating formula was a fully formulated powder containing Eudragit L100-55 (supplied by Degussa) and other pharmaceutically acceptable additives including Yellow#6 lake pigment.

The coating dispersions were prepared using distilled water.

### **Experimental Design Software**

The experimental design and data analysis were conducted using D.o.E. Fusion software supplied by S-Matrix Corp. (Eureka, CA).

### **Experimental Design Parameters**

 Independent Variables: Exhaust Heat (°C) 30-35
 Fluid Delivery Rate (grams/minute) 60-80
 Pan Speed (rpm) 10-15
 % Solids (in Water) 20-30

Number of Experimental Runs: 32 with 4 replicate pairs

### **Dispersion & Coating**

 Dispersion Preparation: Silverson Hi-Shear Mixer: L4RT-A RPM: 10,000 Duration: 10 min/run

- Response Variables
   Tablet Surface Roughness
   Aspirin Dissolution
   Standard Enteric Disintegration
   Stressed Enteric Disintegration
   Coating Process Efficiency
- Coating Equipment/Conditions: 24" O'Hara Labcoat II Spray Guns: 2 Spraying Systems Co (1/8-VAU-SS) Air caps: VA1282125-60-SS Atomization Air: 35 psi (2.4 bar) Pattern Air: 35 psi (2.4 bar) Pan Charge: 14 kg Inlet Air Flow Rate: 250 ft.<sup>3</sup>/minute (7.1m<sup>3</sup>/min)

Tablet Samples Used:

- $\Rightarrow$  325 mg
- $\Rightarrow$  Sub-coat: Opadry II<sup>®</sup>, high performance film coating system, Y-30-18037 (theoretical 2% weight gain)
- Acrylic Coating Weight Gain: theoretical 10% for each run



## **Analytical Methods**

Dissolution:	USP "Delayed-Release Aspirin Tablet" Monograph Q (pass) = 80% dissolved in 90 minutes Results expressed as time required for 80% of the aspirin to dissolve ("DT-80%")
Enteric Disintegration (ET):	Modified USP <701> using 50 rather than 6 tablets pH 1.0 (0.1 N HCl) for 1 hr Pass = no signs of disintegration
Stressed Enteric Disintegration (SET):	Friabilation (50 Tabs, 4 min, 100 revolutions) ET protocol subsequently followed
Surface Roughness:	Determined by MicroPhotonics (Allentown, PA) Average roughness (Ra) determined on five tablets/run
Coating Process Efficiency:	% Efficiency = <u>100*(theoretical - actual weight gain)</u> theoretical weight gain

# RESULTS

Run No.	Exhaust Heat	FD Rate (g/min)	Pan Speed (rpm)	% Solids	DT-80% (min)	Coating Efficiency (%)	SET % Pass	Roughness
1	30	60	15	20	48	75	68	0.91
2 (R1)	33	70	13	25	36	74	64	1.19
3	35	70	13	25	39	74	44	1.2
4	30	80	10	30	46	85	80	1.47
5	35	60	15	20	48	76	76	1.02
6	30	80	15	30	41	82	74	1.21
7 (R2)	35	80	15	25	43	78	38	1.24
8	34	75	14	28	44	76	42	1.17
9	35	80	15	20	42	77	72	1.04
10 (R3)	35	70	15	30	25	73	34	1.32
11	35	80	15	30	32	75	64	1.47
12 (R3)	35	70	15	30	27	69	36	1.31
13	33	70	10	25	28	72	70	1.27
14 (R4)	35	60	10	30	24	65	22	1.43
15	30	80	10	20	42	77	98	1.06
16	35	80	13	30	40	72	54	1.48
17	34	65	11	23	53	64	52	1.18
18	35	80	10	30	41	68	42	1.51
19	30	60	15	30	32	67	18	1.21
20	35	60	15	30	28	64	14	1.29
21 (R4)	35	60	10	30	27	71	26	1.52
22	35	60	10	20	39	75	74	1.18
23	31	65	11	28	40	79	54	1.25
24 (R1)	33	70	13	25	40	78	68	1.2
25	30	70	13	25	41	81	74	1.09
26	35	80	10	20	46	79	92	1.18
27	30	80	15	20	57	90	100	0.92
28	30	60	10	20	38	79	88	1.05
29	34	65	14	23	32	76	50	1.13
30 (R2)	35	80	15	25	41	79	38	1.22
31	34	70	12	30	25	68	74	1.51
32	34	70	12	25	28	70	90	1.28

### Table 1. Replicate Pairs Indicated by R1-R4 Designations




### **General Observations**

- ET performance was insensitive to process conditions.
  - $\circ$  98-100% of the tablets from each run passed
- All tablets were smooth; however, discernible differences in roughness were measured.
- No tackiness was observed in any of the runs.
- The standard deviation for each response variable measurement was < 5% based on replicate pair data analysis.


### Key for Graphs:

- LB = lower bound for selected independent variable
- UB = upper bound for selected independent variable



### Figure 1. Coating Efficiency vs. % Solids and Fluid Delivery Rate

Figure 2. Coating Efficiency vs. Pan Speed and Exhaust Heat SET Performance vs. % Solids and Fluid Delivery Rate





Colorcon





Figure 3. SET Performance vs, Solids and Fluid Delivery Rate

Figure 4. SET Performance vs, Pan Speed and Exhaust Temperature







Figure 6. Fluid Delivery Rate Surface Roughness vs. Pan Speed and Exhaust Heat









Figure 7. Dissolution Time vs. Pan Speed and Exhaust Heat

Figure 8. Dissolution Time vs. % Solids and Fluid Delivery Rate



Colorcon

### Table 2. Coating Parameters Optimized for Coating Efficiency (as calculated by D.o.E. Fusion)

Parameter	Value
% Solids	20
Fluid Delivery Rate (g/min)	80
Pan Speed (rpm)	15
Exhaust Temperature (C)	30
Predicted Coating Efficiency Response = 87%	
Actual Coating Efficiency (Run #27) = 90%	
(Parameters used in Run #27 were identical to those listed	d above.)

Response Range		Run # 27
Response Variable	For All Runs	Result
ET (% Pass)	98-100	100
SET (% Pass)	14-100	100
DT-80% (Minutes)	25-57	57
Roughness (Ra)	0.91-1.52	0.92
Coating Efficiency (%)	64-90	90

## CONCLUSIONS

Coating process efficiency is maximized when:

- % solids and exhaust temperature are minimized
- Fluid delivery rate and pan speed are maximized

The SET pass rate and dissolution time were also maximized under the same conditions, since the amount of polymer applied to each tablet increased as the coating efficiency increased.

Surface roughness was minimized under these conditions due to the reduced incidence of spray drying.

The developmental film coating formula can be applied under high productivity conditions to obtain tablets which pass the USP Delayed Release Aspirin criteria.

Reprint of poster presented at AAPS - Nov 2000.



For more information, contact your Colorcon representative or call:

 North America
 Europe/Middle East/Africa
 Asia

 +1-215-699-7733
 +44-(0)-1322-293000
 +65

Asia Pacific **+65-6438-0318**  Latin America +54-11-4552-1565 © Colorcon, 2009. The information contained in this document is proprietary to Colorcon and may not be used or disseminated inappropriately.

All trademarks, except where noted, are property of BPSI Holdings, LLC.

You can also visit our website at www.colorcon.com