

## Evaluation of the Enteric Performance of Lansoprazole Mini-Tabs Coated in a Perforated Pan

### PURPOSE

There is a growing interest in multi-particulate (MP) modified release (MR) drug delivery systems due to their consistent and reliable in-vivo drug release, with a reduced risk of local irritation along the gastrointestinal tract. <sup>(1, 2)</sup>

Mini-tab technology combines the advantages of MP dosage forms with established manufacturing techniques used in tableting. Additional benefits of mini-tabs include excellent size uniformity, regular shape and a smooth surface, thereby offering an excellent substrate for coating with MR polymeric systems.

The aim of this study was to investigate the potential of successful manufacture and enteric coating of lansoprazole mini-tabs using methacrylic acid co-polymers (Acryl-EZE® aqueous acrylic enteric system, and Eudragit L 30 D-55) in a perforated pan. The produced enteric protection was assessed by both “acid-uptake” test and USP drug release criteria. Lansoprazole is a proton pump inhibitor (PPI), and therefore for the “acid-uptake” test, acid resistance of the enteric coated mini-tabs in 0.1N HCl and at intermediate pH conditions (acetate buffer pH 4.5) was investigated. <sup>(3)</sup>

### METHODS

#### *Formulation and Manufacture of Mini-Tabs*

All ingredients (Table 1) except for the magnesium stearate were screened through a 35 mesh (500 µm) sieve and blended in a Turbula mixer for 10 minutes. Magnesium stearate was then added and the complete formulation was blended for an additional minute.

Mini-tabs were produced by direct compression on a 10-station instrumented Piccola press (Riva, Argentina) using 2 mm standard concave tooling (Notter GmbH, Germany) at 35 rpm and 1.2 kN compression force.

**Table 1. Lansoprazole Mini-Tab Formulation**

<b>Material</b>	<b>% w/w</b>	<b>mg/tablet</b>
Lansoprazole (Jenson Pharm. Services, UK)	12.0	0.98
Lactose (Fast Flo, Kerry Bio-Science, Ireland)	61.4	5.03
Magnesium carbonate (VWR International Ltd., UK)	13.0	1.07
StarCap 1500® co-processed starch excipient (Colorcon)	13.0	1.07
Magnesium stearate (Peter Greven, Holland)	0.5	0.04
Fumed silica (Aerosil 200, Evonik, Germany)	0.1	0.01
Total	100.0	8.20

### Coating of Mini-Tabs

Coating trials were conducted in an Labcoat II-X (O'Hara, Canada) fitted with a custom-made 10" perforated pan with 400 mm mesh insert, a 1 mm nozzle ABC spray gun (Schlick, Germany) and a Perspex shield plate (Colorcon) positioned at the front of the pan to prevent mini-tab loss during coating. The mini-tabs (500 g batch size) were seal-coated with an Opadry® complete film coating system, clear, YS-1-7006 to 5% weight gain (WG) in order to improve mechanical strength of the substrate prior to application of the functional coat. Then the mini-tabs were enteric-coated with either Acryl-EZE, white (93A18597) or L30 D-55 (Evonik, Germany) up to 35% WG.

### Testing of Lansoprazole Mini-Tabs

The mechanical strength of ten uncoated and coated mini-tabs was determined using Schleuniger-4M hardness (Schleuniger, Switzerland) and Copley TA friability (Erweka GmbH, Germany) testers.

Enteric protection of the Acryl-EZE coated mini-tabs was evaluated by an "acid uptake" test. Enteric coated mini-tabs (n=6) were weighed individually and placed in a Copley ZT54 (Copley, UK) disintegration tester with a round 500 µm mesh fitted to the base of the testing tubes. After 2 hours in either 0.1N HCl or pH 4.5 acetate buffer at 37 ± 1°C, the mini-tabs were removed from the vessel, excess surface moisture was eliminated with a paper towel and the tablets were re-weighed. The difference in the respective weights was reported as percent acid uptake.


Dissolution testing of the enteric coated mini-tabs in hard gelatin capsules was conducted in a dissolution bath (SOTAX, Switzerland) according to the USP monograph for "Lansoprazole Delayed-Release Capsules" <sup>(4)</sup>, using Apparatus II (paddles) at 75 rpm. Drug release was measured for a 15 mg lansoprazole dose, i.e. 15 mini-tabs were placed into a size 3 hard gelatin capsule. Dissolution protocol was 1 hour in 0.1N HCl (acid stage) followed by pH 6.8 phosphate buffer with sodium lauryl sulphate (buffer stage) at 37 ± 0.5°C. A dual beam UV/VIS spectrophotometer (Lambda 25, Perkin Elmer Instruments, UK) was used for lansoprazole detection at a wavelength of 306 nm in acid and at 285 nm in buffer. (4) Measurements at each time point were performed in triplicate, and mean and standard deviation (SD) values were calculated.

## RESULTS

Coated mini-tabs exhibited good appearance, showing no visual defects in the coating (Figure 1).

**Figure 1. Appearance of Mini-Tabs Compared to 10 mm Diameter Round Tablets**


[(a) 10 mm uncoated; (b) 10 mm enteric-coated, 10% WG; (c) 2mm, uncoated; (d) 2 mm seal-coated; (e) 2 mm, seal-coated and 30% WG enteric-coated]



The mechanical strengths of the mini-tabs improved significantly after the application of a seal-coat. Breaking force increased from 1.1 kp (uncoated) to 3.3 kp (seal-coated) and friability was reduced from 0.24% to less than 0.01% respectively.


Acid uptake results showed that 30-35% Acryl-EZE enteric coating weight gain was required to obtain sufficient acid resistance, in both 0.1N HCl (pH 1.2) and acetate buffer (pH 4.5) media (Figure 2), and acceptable drug release in USP phosphate buffer (pH 6.8).

**Figure 2. Acid Uptake for Acryl-EZE® Coated Mini-Tabs in 0.1N HCl and Acetate Buffer (n = 6)**



Lansoprazole release profiles from the mini-tabs coated with Acryl-EZE, even at a low weight gain of 15%, was found to comply with the USP specification of <10% drug release after 1 hour in 0.1N HCl and >80% after 1 hour in pH 6.8 phosphate buffer (Figure 3).

**Figure 3. Drug Release from Acryl-EZE® Coated Mini-Tabs**




It was also found that after exposure to the acid media, Acryl-EZE coated mini-tabs were sticking to each other (Figure 4) leading to slower and variable drug release profiles when tested in buffer (Figure 3). These results confirm the findings of Deng et al. (5) who reported similar phenomenon with pellets coated with enteric systems based on methacrylic acid copolymer.

To further investigate the mini-tab sticking issue in the dissolution apparatus, the lansoprazole mini-tabs were enteric coated with Acryl-EZE or Eudragit L 30 D-55. Agglomeration of mini-tabs was observed for both systems after 1 hour in 0.1N HCl, affecting their dissolution profiles in buffer phase (Figure 5). Further work is underway to investigate possible methods of preventing mini-tabs agglomeration in acid.

**Figure 4. Agglomeration of Enteric-Coated Mini-Tabs During Dissolution Testing**



**Figure 5. Drug Release from Mini-Tabs Coated with Eudragit L 30 D-55 or Acryl-EZE®**



## CONCLUSIONS

Lansoprazole mini-tabs (2 mm diameter) were successfully manufactured and enteric coated using a custom-made 10"perforated coating pan, with excellent coating appearance attained. The mechanical strength of the mini-tabs significantly improved upon application of a seal-coat. The enteric coated mini-tabs exhibited low acid uptake in pH 1.2 HCl and pH 4.5 acetate buffer, and consistent drug release in USP phosphate buffer pH 6.8.

Future studies will investigate the mechanism of agglomeration during dissolution testing and potential resolution of this issue related to mini-tabs coated with enteric acrylic systems.

Reprint of poster presented at AAPS 2008.

Authors: Hua Vuong, Marina Levina and Ali Rajabi-Siahboomibi (Colorcon)

## REFERENCES

1. Khosla R., Feely L.C., Davis S.S., Gastrointestinal transit of non-disintegration tablets in fed subjects. *Int. J.Pharm.*, 53 (1989) 107-117.
2. Riis T., Bauer-Brandl A., Wagner T., Kranz H., pHindependent drug release of an extremely poorly solubleweakly acidic drug from multi-particulate extended release formulations. *Eur. J. Pharm/ Biopharm.*, 65 (2007) 78-84.
3. Rohss K., Lind T., Wilder-Smith C. Esomeprazole 40mg provides more effective intragastric acid control than lansoprazole 30 mg, omeprazole 20 mg, pantoprazole
4. 40 mg and rabeprazole 20 mg in patients with gastroesophageal reflux symptoms. *Eur. J. Clin. Pharmacol.*,60 (2004) 531-9.
5. USP 27-NF22 (2004), Lansoprazole Delayed-Release Capsules, pp. 1069-1070.
6. Deng H., Reyes G., Vass S., Rajabi-Siahboomi A.R.,In-vitro dissolution testing of delayed release multi-particulate systems, CRS 2008.

---

The information contained herein, to the best of Colorcon, Inc.'s knowledge is true and accurate. Any recommendations or suggestions of Colorcon, Inc. with regard to the products provided by Colorcon, Inc. are made without warranty, either implied or expressed, because of the variations in methods, conditions and equipment which may be used in commercially processing the products, and no such warranties are made for the suitability of the products for any applications that you may have disclosed. Colorcon, Inc. shall not be liable for loss of profit or for incidental, special or consequential loss or damages.

Colorcon, Inc. makes no warranty, either expressed or implied, that the use of the products provided by Colorcon, Inc., will not infringe any trademark, trade name, copyright, patent or other rights held by any third person or entity when used in the customer's application.

For more information, contact your Colorcon representative or call:

North America  
+1-215-699-7733

Europe/Middle East/Africa  
+44-(0)-1322-293000

Asia Pacific  
+65-6438-0318

Latin America  
+54-11-5556-7700

You can also visit our website at [www.colorcon.com](http://www.colorcon.com)



© BPSI Holdings LLC, 2013.

The information contained in this document is proprietary to Colorcon and may not be used or disseminated inappropriately.

All trademarks, except where noted, are property of BPSI Holdings, LLC.

PR\_Vuong\_mini-tabs\_per\_pan\_ve3\_12\_2013