

Application of Powder Layering Technology and Aqueous Enteric Coating of Lansoprazole 15mg Pellets

OBJECTIVES

Lansoprazole delayed release pellets have been shown to have better absorption properties than a delayed release tablet.¹ However, binder solutions for drug layering and extrusion-spheronization excipients are incompatible with lansoprazole.² In contrast, dry powder layering technology has been reported to provide a more stable manufacturing method for acid-labile drugs.³ The objective of this study was to evaluate the performance of lansoprazole powder layered pellets coated with Acryl-EZE[®], aqueous acrylic enteric system, 93F19255 in various media.

METHODOLOGY

Powder Layering

Lansoprazole was applied to sugar spheres (840-1000 μ m) in a centrifugal fluid bed granulator (Glatt, GPCG-1) using dry powder layering technology. The dusting powders (Table 1) were applied individually while spraying a 3% HPC binder solution (approximately 1340g and 660g, respectively). The second dusting powder was applied to the pellets to enhance drug protection. The process parameters used are listed in Table 2.

Aqueous Enteric Coating

The drug layered pellets were dried (40°C) and screened (1190 μ m) prior to coating in a Strea-1 Aeromatic fluid bed coater (Schlick gun) with Acryl-EZE 93F19255 to a 26% theoretical weight gain. The coating dispersion was prepared under low shear mixing (Table 3), screened (250 μ m), and applied to the pellets using the processing parameters listed in Table 4.

Dissolution Testing

Dissolution testing was performed in a USP apparatus 2 (VanKel VK7000 dissolution bath) at 75rpm, 37.0 \pm 0.5°C. The delayed release dissolution testing (n=6) was performed in 500ml of 0.1N HCl or 0.05M sodium acetate buffer (pH 4.5), followed by phosphate buffer (pH 6.8). Lansoprazole absorbance was measured at 306nm (Agilent 8453 spectrophotometer).

Acid Phase

Six capsules or bead samples (15mg/253mg beads) were placed in 0.1N HCl or acetate buffer (pH 4.5) for 60 minutes. A 25ml aliquot was removed and scanned for the amount of drug released in the acid phase. The specification for the acid phases was not more than 10% drug dissolved after 60 minutes.

Buffer Phase

25ml of 0.1N HCl was removed from each vessel and 425ml of 0.05M phosphate buffer (pH 6.8) was added to obtain 900ml of pH 6.8 media, or all of the 0.05M acetate buffer was removed and 900ml of 0.05M phosphate buffer (pH 6.8) was added to obtain a final pH of 6.8. Sample aliquots were withdrawn from the dissolution vessels at 5, 10, 15, 30, 45, and 60 minutes and analyzed for the amount of drug dissolved. The specification for the buffer phase was not less than 80% drug dissolved after 60 minutes.

Stability Evaluation

Hand-filled capsules (Capsugel hard gelatin, 15mg/253mg beads) were stored (HDPE bottles) at 30°C/60%RH and 40°C/75%RH for 1 month.

Table 1. Lansoprazole Powder Layering Formulation⁴

| Step #1: Dusting Powder | Function | Supplier | % w/w |
|--|---------------------|-------------------------|-------|
| Lansoprazole | Active | Cadila Ahmedabad, India | 13.0 |
| Magnesium Carbonate, Heavy | Stabilizer | EDM Chemical NJ, USA | 13.0 |
| Sucrose | Filler | Domino Sugar NY, USA | 13.0 |
| Corn Starch | Filler | Staley Starch IL, USA | 13.0 |
| L-HPC | Disintegrant Binder | Biddle Sawyer NY, USA | 13.0 |
| Step #: Dusting Powder | | | |
| Sucrose | Filler | Domino Sugar NY, USA | 12.2 |
| Corn Starch | Filler | Staley Starch IL, USA | 10.4 |
| L-HPC | Disintegrant Binder | Biddle Sawyer NY, USA | 10.4 |
| Steps #1 and 2: Binder Solution | | | |
| HPC and Water | Binder | Shin Etsu Tokyo, Japan | 3.0 |
| Total (solids): | - | - | 100 |

Table 2. Powder Layering Parameters

| Parameters | Value |
|-------------------------------------|---------|
| Batch Size (g), 840-100µm nonpareil | 2250 |
| Rotor Speed (rpm) | 200 |
| Binder Spray Rate (g/min) | 20 |
| Powder Addition Rate (g/min) | 15 |
| Inlet Air Temperature (C) | 55 |
| Outlet Air Temperature (C) | 45 |
| Bed Temperature (C) | 45 |
| Atomization Air Pressure (bar) | 1.5 |
| Air Flap (%) | 20 |
| Air Flow (m ³ /hr) | 68 – 80 |
| Total Processing Time (min) | 113 |

Table 3. Dispersion Preparation

| Parameters | Acryl-EZE 93F19255 |
|-------------------------------|--------------------|
| Dispersion Solids Content (%) | 20 |
| Theoretical Weight Gain (%) | 26 |
| Powder (g) | 130 |
| Deionized Water (g) | 520 |
| Total Dispersion (g) | 650 |
| Dispersion Mixing Time (min.) | 30 |

Table 4. Enteric Coating Parameters

| Parameters | Value |
|---|-------|
| Batch Size (g). 1190µm Drug Layered Pellets | 500 |
| Coating Spray Rate (g/min) | 4.5 |
| Inlet Air Temperature (C) | 44 |
| Outlet Air Temperature (C) | 33 |
| Bed Temperature (C) | 32 |
| Atomization Air Pressure (bar) | 1.2 |
| Total Processing Time (min) | 144 |

RESULTS

Lansoprazole powder layered pellets prepared in a centrifugal fluid bed granulator, increased in size from 840-1000µm (nonpareils) to 1190-1410µm. No purple/brown discoloration (sign of degradation) was observed. Pellets appeared spherical, dense, had very few fines, and were therefore suitable for aqueous enteric coating.

Figure 1. Drug Layered and Coated Pellets



Lansoprazole pellets coated with Acryl-EZE 93F19255 (26% theoretical weight gain) exhibited excellent enteric protection in acidic media (pH 1.2 and 4.5) and rapid drug release in phosphate buffer (pH 6.8). The results showed no drug release in 0.1N HCl, followed by 80% release in phosphate buffer (pH 6.8) in 15 minutes or no drug release in acetate buffer (pH 4.5) after 60 minutes followed by 80% release in 20 minutes in phosphate buffer (pH 6.8).

Figure 2. Lansoprazole Dissolution Profiles

Dissolution in phosphate buffer (pH 6.8) after 1 hour in 0.1N HCl or acetate buffer (pH 4.5), n=6, Gelatin Capsule




Table 5. Stability Evaluation

| Lansoprazole Pellet Assay | % Label Claim |
|--------------------------------------|----------------------|
| 1 month 30°C/60%RH with desiccant | 97.2 |
| 1 month 30°C/60%RH without desiccant | 97.2 |
| 1 month 40°C/75%RH with desiccant | 96.0 |
| 1 month 40°C/75%RH without desiccant | 96.2 |

CONCLUSIONS

Dry powder layering technology was used to prepare pellets of an acid-labile drug, lansoprazole (15mg), suitable for aqueous enteric coating. Acryl-EZE 93F19255 provided enteric protection in 0.1N HCl and acetate buffer (pH 4.5) and immediate drug release in phosphate buffer (pH 6.8).

Reprint of poster presented at AAPS – Nov 2006. Authors: Cara J. Young, Kurt A. Fegely, and Ali R. Rajabi-Siahboomi.

REFERENCES

1. Tabata, T., et al., Journal of Biopharmaceutical Sciences, 2 (4), 319-328 (1991)
2. Tabata, T., et al., Drug Development, Industrial Pharmacy, 18 (13), 1437-1447 (1992)
3. Tabata, T., et al., Drug Development, Industrial Pharmacy, 20 (9), 1661-1672 (1994)
4. United States Patent 5,045,321

For more information, contact your Colorcon representative or call:

| | | | |
|------------------------|----------------------------|----------------------|-------------------------|
| North America | Europe/Middle East/Africa | Asia Pacific | Latin America |
| +1-215-699-7733 | +44-(0)-1322-293000 | +65-6438-0318 | +54-11-4552-1565 |

You can also visit our website at www.colorcon.com



© Colorcon, 2009. The information contained in this document is proprietary to Colorcon and may not be used or disseminated inappropriately.

All trademarks, except where noted, are property of BPSI Holdings, LLC.

ads_acryleze_app_powd_lay_v2_07_2009

This document is valid at the time of distribution. Distributed 24-Jan-2021 (UTC)