

# Investigation of the Effect of Ethylcellulose Viscosity Variation Using QbD Samples on Drug Release from Extended Release Multiparticulate

CRS Poster Reprint  
2014

Raxit Y. Mehta<sup>1</sup>, Jason Teckoe<sup>1</sup>, Cody Schoener<sup>2</sup>, Shari Workentine<sup>2</sup> and Ali R. Rajabi-Siahboomi<sup>1</sup>

<sup>1</sup>Colorcon Inc., Harleysville, PA 19438, USA  
www.colorcon.com

## Abstract Summary

Extended release multiparticulates (MP) of acetaminophen (APAP) were prepared using ETHOCEL™ Premium Ethylcellulose Polymers (International Flavors and Fragrances Inc., USA) as a rate controlling barrier membrane coating polymer. Quality by Design (QbD) samples of ETHOCEL Std. 10, 20 and 100 Premium grades were organically coated onto drug layered APAP MP up to a 15% weight gain (WG). The objective of the study was to investigate the effect of viscosity variation, within the manufacturer's specification of ETHOCEL Std. Premium grades, on in vitro drug release from APAP MP.

## Introduction

Due to its excellent barrier properties, safety and global regulatory acceptance, ethylcellulose is a commonly used polymer for extended release (ER) MP formulations.<sup>1</sup> The rate of drug release from ER MP is essentially controlled by the properties of the barrier membrane coating.<sup>2,3</sup> The critical material attributes of ethylcellulose that may influence drug release were considered to be viscosity and ethoxyl substitution. The aim of this study was to evaluate the effect of viscosity variation for ETHOCEL Std. 10, 20 and 100 Premium grades on the drug release from APAP ER MP. Samples of ethylcellulose coated APAP MP were tested at 5% and 15% WG to assess drug release variability.

## Experimental Methods

### Drug Layering

The composition of the drug layered MP is displayed in Table 1. Uncoated sugar spheres (SUGLETS® #18-20 Colorcon Inc., USA) were coated in an Oyster Huttlin Unilab fluid bed (Huttlin GmbH, Germany) using semi-fine grade of APAP (Mallinckrodt, USA) and HPMC based Opadry® (Colorcon Inc., USA) as a binder. Drug to binder ratio of 50:50 w/w was maintained. The drug layered MP were screened to remove agglomerates and fines before application of the barrier membrane coating.

Table 1. Composition of Drug Layer

,QJUHGLHQWV	6XSSOLHU	% Z/Z
\$FHWDP LQR SKHQ 863	&RYLGLHQ, 86\$	7
6HPL-ILQH 3RZGHU		
850/ 1000 6XJOHWV (6XJDU6SKHUHV1) 18/20	&RORUFRQ,86\$	86
2SDGU\	&RORURQ, 86\$	7
7RWDO		100

Table 2. ETHOCEL Sample Information

9LVFRVLW\ *UDGH	9LVFRVLW\ F3	4E' 6DP SOH F3	&RDWLQJ 6ROLGV % (Z/Z)
(7+2&/ 6WG. 10 3UHPLXP	9-11	9, 10, 11	7
(7+2&/ 6WG. 20 3UHPLXP	18-22	18, 20, 22	5
(7+2&/ 6WG. 100 3UHPLXP	90-110	90, 100, 110	3

Table 3. ETHOCEL Coating Process Parameters

3URFHVV 3DUDPHWHU	9DOXH
%DWFK VLJH (J)	750
,QOHW WHP SHUDWXUH (R&)	38 ± 42
3URGXFW WHP SHUDWXUH (R&)	30 ± 32
2XWOHW WHP SHUDWXUH (R&)	29 ± 31
\$WRPLJLQJ DLU (EDU) / (SVL)	1.3 / 18.8
\$LU YROXPH (P <sup>3</sup> /KU) / (FIP)	45 - 50 / 26.5 - 29.5
)OXLG GHOLYHU\ UDWH (J/PLQ)	5 ± 7
&RDWLQJ VROXWLRQ YLFRVLWV (F3)	70-85

### Application of Ethylcellulose

The ETHOCEL QbD samples used in the study are shown in Table 2. These samples cover the viscosity specification range for each ETHOCEL Std. Premium grade. Dibutyl sebacate (Vertellus, USA) was used as a plasticizer at a 9:1 w/w ratio of polymer to plasticizer. The organic coating solutions of ethylcellulose were prepared using isopropyl alcohol and purified water (90:10 w/w) as the solvent. Organic coating trials were carried out using a Glatt GPCG-2 (Glatt Air Techniques Inc., USA) fluid bed. The process parameters for applying the ethylcellulose coatings are listed in Table 3.

### Dissolution Studies

In vitro dissolution studies were carried out using USP Apparatus I (baskets) at 100 rpm in 1000 ml of purified water. Drug release was determined spectrophotometrically at a wavelength of 243 nm. Drug release data for all MP were compared using similarity factor ( $f_2$ ) analysis.

## Results and Discussion

The drug release profiles for MP coated with ETHOCEL Std. 10 Premium QbD samples are shown in Figures 1 and 2 at 5% and 15% WG, respectively. The ( $f_2$ ) similarity factor values were  $f_2 > 75$  at a 5% WG and  $f_2 > 83$  at 15% WG. This indicates there is minimal variability of APAP drug release due to variation of viscosity within the ETHOCEL Std. 10cP Premium grade specification.

Figure 1: Drug Release from ETHOCEL Std. 10 Premium QbD Samples at 5% WG






Figure 2: Drug Release from ETHOCEL Std. 10 Premium QbD Samples at 15% WG




The drug release for ETHOCEL Std. 20 Premium QbD samples are shown in Figures 3 and 4, while the MPs coated with ETHOCEL Std. 100 Premium QbD samples are shown in Figures 5 and 6.

**Figure 3: Drug Release from ETHOCEL Std. 20 Premium QbD Samples at 5% WG**



**Figure 4: Drug Release from ETHOCEL Std. 20 Premium QbD Samples at 15% WG**



The drug release  $f_2$  values for 5% WG and 15% WG ETHOCEL Std. 20 Premium samples were  $f_2 > 74$  and  $f_2 > 64$  respectively, while the equivalent coating weight gain of the ETHOCEL Std. 100 Premium samples were  $f_2 > 77$  and  $f_2 > 70$  respectively. This indicates there is minimal variability of APAP drug release due to variation of viscosity within the ETHOCEL grade specifications for both ETHOCEL Std. 20 Premium and ETHOCEL Std. 100 Premium.

In all cases, increased weight gain of the ethylcellulose barrier membrane coating resulted in slower drug release, while use of higher viscosity grades resulted in longer initial lag times.

**Figure 5: Drug Release from ETHOCEL Std. 100 Premium QbD Samples at 5% WG**





Figure 6: Drug Release from ETHOCEL Std. 100 Premium QbD Samples at 15% WG



### Conclusions

This study shows that viscosity variation, within the manufacturer’s specifications for ETHOCEL Std. 10, 20 and 100 Premium grades, has minimal impact on drug release for APAP ER MP at both low (5%) and high (15%) weight gains. These results highlight the consistency of the ETHOCEL product and utility of ETHOCEL QbD samples as a means to develop robust ER MP dosage forms.

### References

1. ETHOCEL Product Information Brochure, accessed on January 17, 2014  
<http://www.colorcon.com/literature/marketing/mr/Extended%20Release/ETHOCEL/Chinese/730%20Colorcon%20Ethocel2%20w-new%20logo-09.03.09.pdf>
2. Missaghi, S. et al., Controlled Release Society Annual Meeting, New York, NY, July 2008.
3. Yu, L. X., Pharmaceutical Quality by Design: product and process development, understanding, and control. Pharmaceutical Research 2008;25(4): 781-791.

The information contained herein, to the best of Colorcon, Inc.’s knowledge is true and accurate. Any recommendations or suggestions of Colorcon, Inc. with regard to the products provided by Colorcon, Inc. are made without warranty, either implied or expressed, because of the variations in methods, conditions and equipment which may be used in commercially processing the products, and no such warranties are made for the suitability of the products for any applications that you may have disclosed. Colorcon, Inc. shall not be liable for loss of profit or for incidental, special or consequential loss or damages.

Colorcon, Inc. makes no warranty, either expressed or implied, that the use of the products provided by Colorcon, Inc., will not infringe any trademark, trade name, copyright, patent or other rights held by any third person or entity when used in the customer’s application.

For more information, contact your Colorcon representative or call:

North America  
+1-215-699-7733

Europe/Middle East/Africa  
+44-(0)-1322-293000

Asia Pacific  
+65-6438-0318

Latin America  
+54-1-5556-7700

All trademarks, except where noted, are property of BPSI Holdings LLC. The information contained in this document is proprietary to Colorcon, Inc. and may not be used or disseminated inappropriately.

ETHOCEL™ is a trademark of International Flavors and Fragrances Inc. or its affiliates. © 2021 IFF. All rights reserved

©BPSI Holdings LLC 2014