

Coating Parameters – Pigmented PVA-Based Formulations

Nutrafinish[®], dietary supplement coatings are one-step film coating products available in clear, white and pigmented formulations. Developed specifically for food supplements, they provide excellent film forming capabilities including logo definition, high tensile strength and excellent adhesion properties.

The coating parameters listed in the table below are based on Colorcon data for PVA-based pigmented formulations. Individual product and machine functions should be considered and conditions altered as required.

		Fully Perforated Side Vented Coating Pan				Conventional Pan
Pan Diameter	cm	38.1	60.96	122	152.4	40.64
	in	15	24	48	60	16
			Solution Prepara	tion		
Solvent		Purified Water	Purified Water	Purified Water	Purified Water	Purified Water
Solids Content	% w/w	20	20	20	20	20
			Run Objective	S		
Theoretical Weight Gain	% w/w	3	3	3	3	3
Tablet Charge [†]	kg⁺	1.5 - 3.0	12 - 15	120 - 160	320 - 350	3
			Gun Configurat	ion		
Example Spray Gun	type	Schlick 931/7-1 S22	Spraying Systems 1/8 VAU	Schlick 931/7-1 S35 (ABC)	Manesty	Schlick 970/7-7 S75
Number of Guns	guns	1	1 - 2	3 – 4	5 - 6	1
Nozzle	mm	1.2	1.2	1.2	1.2	1.2
Diameter	in	0.05	0.05	0.05	0.05	0.05
Gun to Bed	cm	10	17	25	25	15
Distance	in	4	7	10	10	6
Gun to Gun	cm	N/A	10	15	15	N/A
Distance	in	N/A	4	6	6	N/A
Atomizing Air	bar	1.4	1.7	2.5	2.5	1.4
Pressure	psi	20	25	36	36	20
Pattern Air	bar	1.4	1.7 25	2.5	2.5	1.4
Pressure	psi	20		36	36	20
			Process Parame			
Drying Air	cu.m/hr	250 - 290	424 - 470	2300 - 3060	5440 - 6630	140
Volume	cu.ft/min	147 - 175	250 - 275	1350 – 1800	3200 - 3900	80
Pan Speed*	rpm*	18	14	8	6	20
Pre-warm	°C	45 - 50	45 - 50	45 - 50	45 - 50	45 - 50
Tablet Bed	°F	113 - 122	113 - 122	113 - 122	113 - 122	113 - 122
Inlet Air	°C	60	65	65	65	65
Temperature	°F	140	149	149	149	149
Tablet Bed Temperature	°C	45 - 48	45 - 48	45 - 48	45 - 48	45 - 48
	°F	113 - 118	113 - 118	113 - 118	113 - 118	113 - 118
Spray Rate*	g/min	12 - 20	35 - 50	250 - 350	480 - 680	12

^{*}Pan speed should be set so that even mixing of the bed occurs. This will be influenced by tablet size, shape, friability and baffle type in the pan.

type in the pan.

†Tablet charge is dependent on tablet size and shape.

*Spray rate numbers are estimated.



Nutrafinish®

The information contained herein, to the best of Colorcon, Inc.'s knowledge is true and accurate. Any recommendations or suggestions of Colorcon, Inc. with regard to the products provided by Colorcon, Inc. are made without warranty, either implied or expressed, because of the variations in methods, conditions and equipment which may be used in commercially processing the products, and no such warranties are made for the suitability of the products for any applications that you may have disclosed. Colorcon, Inc. shall not be liable for loss of profit or for incidental, special or consequential loss or damages.

Colorcon, Inc. makes no warranty, either expressed or implied, that the use of the products provided by Colorcon, Inc., will not infringe any trademark, trade name, copyright, patent or other rights held by any third person or entity when used in the customer's copylication.

For more information, contact your Colorcon representative or call:

North America Europe/Middle East/Africa +1-215-699-7733 +44 (0)-1322-293000

Latin America India China +54-11-5556-7700 +91-832-6727373 +86-21-61982300



© BPSI Holding, 2017. The information contained in this document is proprietary to Colorcon, Inc. and may not be used or disseminated inappropriately.

All trademarks, except where noted, are property of BPSI Holdings LLC.

pi_nutrafinish_coat_par_pva_pigm_v3_03_2018



You can also visit our website at www.colorcon.com