

Simplified Formulation of Dietary Supplements Using High Functionality Clean Label Excipients

Authors: Sarah S. Aboelela, Ankit Soni*, Robert Theisen, Amelia Thomas, Rafael Gonzalez, Manish Ghimire, David Ferrizzi

Colorcon, Inc. Harleysville, PA 19438, USA | AAPS Poster Reprint 2023

Introduction

Formulation development of dietary supplements (DS) into robust solid dosage forms can be complex and time consuming. In the nutraceutical sector, consumers desire clean label products containing natural (all plant-based ingredients), easily recognizable ingredients¹. To meet this demand, there is need to develop nutraceutical excipients that reduce time and effort for DS product development.

The purpose of this study was to develop clean label, highly functional excipients for manufacturing of acceptable DS tablets.

- The novel nutraceutical excipients include a Nutracore™ filler excipient blend and Nutracore™ lubricant blend (with both lubricant and glidant properties).
- Three model dietary supplement ingredients were used for tablet formulations: ashwagandha powder, ashwagandha granules and garlic powder.

Objectives:

- To evaluate applicability of clean label Nutracore blends with a range of DS for the powder properties and their tablet performance.
- Competitive evaluation of prepared tablets against two commercially available DS tablets containing either synthetic pharmaceutical excipients or clean label excipients.

Methods

Powder Testing: Properties of blends were evaluated in terms of appearance, bulk density, and particle size distribution (PSD). PSD and flowability were evaluated using Malvern Mastersizer Sirocco 2000 and Flodex™. The three model DS powders were characterized for powder properties and tabletability with Nutracore filler and lubricant blends.

Powder blends were formulated as shown in Table 1 and made using 2-step mixing where the DS and the Nutracore filler blend were mixed for 3 minutes in a V-blender, followed by the lubricant blend and mixing for additional 2 minutes. Tablets were directly compressed into 450 mg, 10 mm round, biconvex tablets using a Piccola tablet press operated at 50 rpm turret speed and 25 kN compression force.

Table 1: Formulation Composition for the Clean Label Tablet Formulations

Excipient	%w/w in formulation		
Ashwagandha Powder	52.0	-	-
Ashwagandha Granules	-	52.0	-
Garlic Powder	-	-	52.0
Nutraceutical Filler Blend	45.5	45.5	45.5
Nutraceutical Lubricant Blend)	2.5	2.5	2.5
Total	100	100	100

Tablet Evaluation: Prepared tablets were evaluated for appearance, tensile strength (USP <1217>) and disintegration time (DT) (USP <701>).

Evaluation vs. Marketed Product: The Nutracore tablet properties (tensile strength and disintegration time) of were compared to two marketed dietary supplement (DS) products. The garlic containing product is an enteric coated tablet. For better comparison of disintegration time, garlic tablet prepared using Nutracore blends was enteric coated using Nutratric® coating at 3% weight gain (WG).

Results

- Nutracore lubricant blend had a light-brown appearance and free flowing powdery form.
- The filler blend was an off-white, free flowing powder with Flodex values of 7-9 mm and compressibility index of 24.6 %.
- The powder properties of pure clean label blends and the three DS are given in Table 2.

Table 2: Powder Properties of the Nutracore Excipient Blends and DS

Material	Bulk Density (g/mL)	Particle Size (µm)		
		Dv10	Dv50	Dv90
Nutraceutical Filler	0.50	31.5	97.5	202.7
Nutraceutical Lubricant	0.26	2.3	14.9	48.2
Garlic powder	0.52	7.4	35.7	123.6
Ashwagandha powder	0.53	7.5	26.1	62.0
Ashwagandha granules	0.37	105.8	246.0	498.3

Figure 1: Appearance of Clean Label Tablets vs. Marketed Products



Marketed Product Evaluation:

Ashwagandha Dietary Supplement Case Study

- The composition and tablet performance of the Nutracore tablets and the marketed products of ashwagandha are summarized in Table 3.
- The Nutracore formulation prepared with minimal ingredients was compared with the marketed tablets. The prepared Nutracore tablets had shiny, uniform-colored surfaces, without any visual defects (Figure 1).
- The Nutracore tablets (uncoated) prepared with either powder or granular ashwagandha as a DS were significantly superior in terms of tensile strength and DT compared to the marketed products as detailed in Table 3.
- The tensile strength of tablets prepared using Nutracore excipients with powdered ashwagandha was >1.6 MPa and DT of <25 minutes which was slightly higher compared to tablets prepared with granular ashwagandha.

Table 3: Composition of Ashwagandha Tablets and Performance: Evaluation vs Marketed Product

Parameters	Clean Label Ashwagandha Powder Tablets	Clean Label Ashwagandha Granules Tablets	Marketed Ashwagandha Tablets
Composition	Ashwagandha powder, Nutracore Filler and Lubricant	Ashwagandha granules, Nutracore Filler, and Lubricant	KSM-66 Ashwagandha® organic extract, organic dextrose, organic guar gum, organic rice extract, organic rice hulls, organic gum arabic, organic sunflower oil, organic clear coating
Tablet properties	450 mg, round biconvex uncoated tablets	450 mg, round biconvex uncoated tablets	850 mg, round biconvex coated tablets
Performance			
Tensile strength (MPa)	2.62 ± 0.18	1.94 ± 0.15	0.65 ± 0.09
DT (min)	23.7	16.1	58.1

Garlic Dietary Supplement Case Study

- The composition and the tablet performance of the Nutracore tablets and the marketed products of garlic is summarized in Table 4.
- Despite the garlic tablets with Nutracore ingredients being uncoated, they had higher tensile strength compared to the marketed tablets containing synthetic excipients as shown in Table 4.
- The disintegration time for both marketed tablets and tablets prepared with Nutracore excipients (coated with Nutratric a nutritional enteric coating) was comparable.

Table 4: Composition of Garlic Tablets and Performance: Evaluation vs. Marketed Product

Parameters	Clean Label Garlic Tablets	Marketed Garlic Tablets
Composition	Garlic powder, Nutracore Filler, Nutracore Lubricant	Garlic, Dicalcium phosphate, Vegetable cellulose, Calcium carbonate, Ethyl cellulose, Calcium silicate, Vegetable magnesium stearate, Titanium dioxide (color), Natural palm leaf glaze, Enteric coating (sodium alginate and stearic acid)
Tablet properties	450 mg, round biconvex enteric coated with Nutratric at 3% WG	~740 mg, enteric coated, round biconvex
Performance		
Tensile strength (MPa)	2.61 ± 0.10 (Uncoated)	1.75 ± 0.12
DT (min)	46.4	52

Conclusions

Nutraceuticals filler and lubricant blends provided simpler, clean label formulations that resulted in robust tablets for ashwagandha and garlic dietary supplements.

Tablet formulations containing label friendly excipients performed significantly better in terms of mechanical strength and disintegration time compared to existing marketed products.

Nutraceuticals filler and lubricant blends are beneficial for nutraceutical applications by reducing time and cost for development.

References

1. Sarah S. Aboelela, Robert Theisen, Manish Ghimire, Ali Rajabi-Siahboomi. The Effect of a Clean-label, Plant-based Natural Glidant on the Flow Properties of Multiple Model Powders compared to a Synthetic Glidant, AAPS 2022.

The information contained herein, to the best of Colorcon, Inc.'s knowledge is true and accurate. Any recommendations or suggestions of Colorcon, Inc. with regard to the products provided by Colorcon, Inc. are made without warranty, either implied or expressed, because of the variations in methods, conditions and equipment which may be used in commercially processing the products, and no such warranties are made for the suitability of the products for any applications that you may have disclosed. Colorcon, Inc. shall not be liable for loss of profit or for incidental, special or consequential loss or damages.

Colorcon, Inc. makes no warranty, either expressed or implied, that the use of the products provided by Colorcon, Inc., will not infringe any trademark, trade name, copyright, patent or other rights held by any third person or entity when used in the customer's application.

Colorcon is a global company located in North America, Europe, Middle East, Africa, Latin America, India, and China.

For more information website at www.colorcon.com