

Acryl-EZE<sup>®</sup>MP is a fully formulated, dry enteric acrylic coating system dispersible in water, for the application of an enteric film coating to multi-particulate solid dosage forms. Combining the benefits of a fully formulated coating system with a globally accepted enteric polymer (EUDRAGIT L100-55\*), Acryl-EZE MP is readily dispersible in water for easy application.

The enteric coating provides consistent, reproducible delayed release profiles.

## Key Characteristics

- Acrylic polymer for proven enteric performance.
- Ready formulated powder.
- Easy to dispense and disperse - only 20 minutes preparation time.

## Applications

- Multi-particulates, such as beads, granules and drug crystals.

## Packaging

- 25kg carton with double polyethylene bag liners and desiccant between bags.

## Recommended Storage Conditions

- Store product below 30°C and less than 75% relative humidity (RH).
- Re-close polyethylene bags (with desiccant) securely before storing.

## Re-evaluation Period (Warranty)

- 12 months from date of manufacture.

\* Methacrylic acid copolymer type C.

For more information, contact your Colorcon representative or call:

North America  
**+1-215-699-7733**

Europe/Middle East/Africa  
**+44-(0)-1322-293000**

Asia Pacific  
**+65-6438-0318**

Latin America  
**+54-11-4552-1565**

You can also visit our website at [www.colorcon.com](http://www.colorcon.com)



© Colorcon, 2009. The information contained in this document is proprietary to Colorcon and may not be used or disseminated inappropriately.

All trademarks, except where noted, are property of BPSI Holdings, LLC.